




**PHARSMART**

SCIENCE & BEAUTY MADE IN ITALY

**Product Certificate KN95 Respirator Masks**

<b>Product Name</b>	<b>KN95 Respirator Masks</b> 
<b>Specification</b>	WCL-0075
<b>Size</b>	13,5 cm x 15,5 cm x 13 cm
<b>Manufacturer:</b>	Jiangsu Weichuangli New Materials Co., Ltd
<b>Inspector</b>	Inspector 002
<b>Composition</b>	Main mask, nose clip and ear loop. The main mask is made of 4-layer filters, non-sterile and disposable
<b>Material</b>	PP nonwoven, Melt-blown fabric
<b>Application</b>	Effectively filter PM2.5, bacteria, virus, dust and smog
<b>Operative Standard</b>	GB2626-2006 Respiratory Protective Equipment-Non-powered air-purifying particle respirator
<b>Restriction</b>	None
<b>Production Date/ Batch No</b>	See the box. Validity 3 years
<b>Storage</b>	Please store the item in a cool and draughtly place away from sunlight, with relatively humidity less than 50% and no corrosive gas

**PHARSMART Srl a Socio Unico**

Sede legale e operativa : Corso Vittorio Emanuele II, 30 - 20122 Milano ITALIA

P. IVA/VAT No. : 09594770969 - REA No. MI-2101110

TEL: + 39 02.36.50.71.98 – E-mail: [info@pharsmart.com](mailto:info@pharsmart.com)



**This product has been tested and proved to comply with the required quality standards, and therefore it has been approved to leave the factory.**

**How to wear a mask:**

- 1. Stretch a mask and check the integrity of its front surface.**
- 2. Position the mask under your chin with the nose-piece up.**
- 3. Extend the mask to cover your mouth and nose, loop the straps over your ears and adjust them to a comfortable state.**
- 4. Use forefingers and middle fingers of both hands to mold the nosepiece to the shape of your nose by pushing downward and inward.**
- 5. Place both hands completely over the mask and perform positive and negative pressure tests to keep your mask from leaking.**

**Please check the package and masks first, if damaged, do not use.**

**Do not wash. Keep dry. If the item is humid, smelly or difficult to breath in, please replace it immediately.**

**The item is disposable and repeated use is not recommended.**

**PHARSMART Srl a Socio Unico**

**Sede legale e operativa : Corso Vittorio Emanuele II, 30 - 20122 Milano ITALIA**

**P. IVA/VAT No. : 09594770969 - REA No. MI-2101110**

**TEL: + 39 02.36.50.71.98 – E-mail: [info@pharsmart.com](mailto:info@pharsmart.com)**